

Rootcare Instruments




Rootcare Instruments

P.o.box: 2406,
Commissioner Road,
Muhammad Pura Sialkot 51310
Pakistan.
Cell: +92 322 5590098
Email: sales@rootcareinstruments.com
Skype: rootcare.ins
Web : www.rootcareinstruments.com

Certification: CE Marking ISO 13485 : 2016



www.rootcareinstruments.com




Lab Instruments



786-1962
Full Arch
Chrome Plated



786-1963
Denture
Chrome Plated



786-1964
Chrome Plated



786-1965
Crowned Bridge
Chrome Plated



786-1966
Full Denture
Chrome Plated



786-1967
Full Denture
Chrome Plated



786-1968
With Magnet




786-1969
With Magnet


Special Features
Lightweight anatomic articulator made of aluminum die-casting. It can be used following the geometrical laws of medium values (Bonwill triangle). Fixing according to centric. Automatic return to the hinge axis with springs. Perfect stability also in reversed position. Easy access and visibility to the lingual surfaces. Interchangeable self centering aluminum mounting plates or magnetic plates. Pure translation movement (key and instruction included). Anti-shock hinge rod in case articulator falls down. Bigger size lock screws for mounting plates.




786-1970



786-1971







786-1988
Zahle



786-1989
Zahle



786-1990
Zahle



786-1991
Hylin



786-1992
Hylin-Ceramica




786-1993
Hylin-Sacred



786-1994
Hylin-Le Cron



786-1995



786-1996
Le Cron

786-1997
Le Cron

786-1998
Vehe

786-1999
Beale

786-2000
Roach


786-2001
Le Cron

786-2002
Le Cron

786-2003
Gritman



Wax and Modelling Carvers



786-2004
Fig.1
P.K. Thomas

786-2005
Fig.2
P.K. Thomas


786-2006
Fig.3
P.K. Thomas

786-2007
Fig.4
P.K. Thomas

786-2008
Fig.5
P.K. Thomas



786-2009
Fig.1
Fahnenstock
cm. 12,5



786-2010
Fig.1
Lessmann
cm. 12,5



786-2011
Fig.1
Con Scodellino
cm. 12,5



786-2012
Fig.2
Fahnenstock
cm. 17



786-2013
Fig.2
Lessmann
cm. 17



786-2014
Fig.2
Con Scodellino
cm. 17




786-2015
cm. 17



786-2016 cm.17
Gritman



786-2017
cm. 19.5



786-2018
cm. 21.5




786-2019
cm. 20.5



786-2020
cm. 19.5



786-2021
cm. 20



786-2022
cm. 22.5

Saw, Blades, Finger Protector
Sterilizing Forceps



786-2023
Saw for Medels

VTS



786-2024/BL-1
0.12x5 mm



786-2025/BL-2
0.25x5 mm



786-2026/BL-3
0.25x3 mm



786-2027


Langenbeck
Finger protectors



786-2028




786-2029



786-2030
Rogge
23 cm



786-2031/1
20 cm



786-2032/2
30 cm