

Rootcare Instruments




Rootcare Instruments


P.o.box: 2406,
Commissioner Road,
Muhammad Pura Sialkot 51310
Pakistan.
Cell: +92 322 5590098
Email: sales@rootcareinstruments.com
Skype: rootcare.ins
Web : www.rootcareinstruments.com

Certification: CE Marking ISO 13485 : 2016

www.rootcareinstruments.com



ORTHODONTICS




786-2033
Waldsachs
15 cm



786-2034
Universal
16 cm




786-2035
Nano
13.5 cm




786-2036
Adener
14.5 cm



786-2037
Modell marburg
Serrated Jaws
14 cm




786-2038
Modell marburg
Smooth jaws
14 cm




786-2039
Modell marburg
Carbide inserts jaws
14 cm




786-2040
Goslee
One Serrated Jaws
14 cm




786-2041
Adorer
12 cm



786-2042
Adorer- Heavy
12 cm




786-2043
Tweed
14.5 cm




786-2044
O'Brien
13.5 cm



786-2045
Wire bending pliers
12.5 cm



786-2046
Adams
12.5 cm



786-2048
Young
14 cm




786-2047
Freyer
12.5 cm



786-2049
Clasp forming
14 cm




786-2050
 Peeso
 13,5 cm



786-2051
 Peeso
 13,5 cm




786-2052
 14 cm



Round nose pliers


786-2053
 Langenbeck
 14 cm



786-2055
 14 cm



786-2056
 13 cm




786-2057
 Angle
 12,5 cm

For wire up to




786-2054
 Andresen
 13 cm



786-2058
 13,5 cm



786-2059
 Schwarz
 13,5 cm



786-2060
 Fischer
 13 cm



786-2061
14 cm



786-2062
14 cm



786-2063
14 cm



786-2064
14,5 cm



786-2065
14,5 cm



786-2066
14 cm




786-2067
14 cm




786-2068
14 cm



786-2069
Serrated Jaws
12 cm



786-2070
How
14,5 cm



786-2071
How-angled
14,5 cm



786-2072
Arrow clasp forming
12 cm



786-2073
Arrow clasp forming
12 cm




786-2074
Arrow clasp bending
12 cm

786-2075

How Plier

Firmly hold the bracket tie wing for seating of the band.
Carbide serrated tips for non slip grip. 49.50



786-2076

De Bonding Plier


Bracket / Adhesive remover with replaceable tool steel tip. Wedge between both edges of the base of the bracket and tooth surface. Easily lifting off the bracket or Bonding material. 49.50



786-2077

Jara Back Plier


A Plier for precise bending of wires. Precision grooves ensure accurate bending and closing of loops. Serrated tip for up to .020" wires. 49.50



786-2078

Posterior Band Remover Long


The padded tip rest on the Occlusal of the molar and the carbide tipped side removes the bracket or Bonding material. 49.50



786-2079

Hard Wire Cutter


Cuts hard wire up to .020" carbide insert for longer life. 49.50



786-2080

Posterior Band Remover Small


49.50



786-2081

Pin & Ligature Cutter


Standard size lips designed to cut soft wire up to .051". 49.50



786-2082

Tweed Loop Forming Plier

3 Size carbide tip allows to make 3 size loops. Can hold wire firmly upto .022" x .028". The other side serrated tip gives none slip grip. 49.50




All Pliers Also Available In Box Joint




786-2083
Bird Beak Plier

1/2" Long tip carbide inserted round tip and flat ground tip gives non slip grip for wire up to .030". 49.50




786-2084
Angled Direct Bond Bracket Remover

60° Angled tip allows for easy access to difficult anterior and posterior areas. 49.50




786-2085
Safety hold, small Distal end Cutter

Narrow and light. 49.50




786-2086
Weingart Utility Plier

Matching deep serrated pads holds wire firmly at a convex working angle. Tapered narrow lip fits easily between brackets. 49.50




786-2087
Safety hold, Distal end Cutter

Carbide tips for longer life and sharpness cuts and hold wire up to .022" x .028". 49.50




786-2088
Three Prong Clasps Adjusting Plier

Rounded jaws close and adjust clasps without niking wire upto .030". 49.50




786-2089
Mini Pin & Ligature Cutter

Narrow tips cuts pins, Ligatures and elastic smoothly and soft wire up to .015". 49.50



786-2090
Light Wire Plier

Wire forming Plier upto .061" the longer tapered beaks make it easier to bend small diameter loops. 3Grooves on square tip. 49.50



All Pliers Also Available In Box Joint



Rectangular Arch plier
13 cm / 5"
786-2091 Box Joint
786-2092 Screw Joint
786-2093 L-Key Joint
Short jaws, Tweed style, for use with square or rectangular wires. Blades parallel at 0.02", 0.51mm opening.



Weingart style plier
13 cm / 5"
786-2094 Box Joint
786-2095 Screw Joint
786-2096 L-Key Joint
Classical design, rounded tips. Useful for removal of ceramics.



Arch contouring plier
13 cm / 5"
786-2097 Box Joint
786-2098 Screw Joint
786-2099 L-Key Joint
For easily forming arches.



Arch forming & contouring plier
13 cm / 5"
786-2100 Box Joint
786-2101 Screw Joint
786-2102 L-Key Joint
With 3 grooves. For wire 0.016", 0.018" & 0.022".





Arch forming plier
13 cm / 5"
786-2103 Box Joint
786-2104 Screw Joint
786-2105 L-Key Joint
Without grooves




Micro Distal end cutter
13 cm / 5"
786-2106 Box Joint
786-2107 Screw Joint
786-2108 L-Key Joint
With safety hold, cuts arch wires. Cuts rectangular wire up to 0.56x0.71mm, round wire up to 0.51mm and twisted wire up to 0.41mm



Hard wire cutter
13 cm / 5"
786-2109 Box Joint
786-2110 Screw Joint
786-2111 L-Key Joint
For cutting orthodontic wires for fixed techniques any size. Cuts rectangular wire up to 0.56x0.71mm and round wire up to 0.51mm



Distal wire cutter
15 cm / 5"
786-2118 Box Joint
786-2119 Screw Joint
786-2120 L-Key Joint



Lindblom
13 cm / 5"
786-2121 Box Joint
786-2122 Screw Joint
786-2123 L-Key Joint



45° Ligature cutter
13 cm / 5"
786-2112 Box Joint
786-2113 Screw Joint
786-2114 L-Key Joint
45° angled for easy access to posterior areas. Cuts ligature wire up to 0.36mm



Miniature ligature cutter
13 cm / 5"
786-2115 Box Joint
786-2116 Screw Joint
786-2117 L-Key Joint
For soft wires. ligature wire up to 0.014" (0.36mm)



Ligature instruments




786-2124
Band seating
&
removing instruments

786-2125

786-2126


786-2127
Elastic Instruments



786-2128
Marshon Band Pusher
14 cm




786-2129
Band Setter with
tin insert
15 cm



786-2130
13.5 cm

786-2131
14 cm



786-2132
Separating pliers for separators
15.5 cm

ORTHODONTIC INSTRUMENTS



786-2133
 Mathieu pliers
 14 cm




786-2134
 Elastic ligature placing pliers
 14 cm




786-2135
 Mathieu ligating pliers
 14 cm



786-2136
 Mathieu plier
 with finger ring
 14 cm



786-2137
 Mosquito plier



786-2139
 Bracket placing
 tweezers



786-2140
 Bracket placing
 tweezers



786-2138
 Bracket positioning gauge



786-2141
 Ligature tying plier




786-2142
 Bracket holding
 tweezers



786-2143
 Band seating
 instrument
 Nylon



786-2144
 Mosquito forceps



Hook

